

Package: tidyprint (via r-universe)

May 30, 2026

Title Custom Print Methods for SummarizedExperiment

Version 1.1.0

Depends R (>= 4.4.0)

Description Provides customized print methods for 'SummarizedExperiment' objects to enhance readability and usability within a tidy workflow. It offers consistent, tidyverse-aligned console displays, including alternative tibble abstractions for large genomic data to improve discoverability and interpretation. The package also includes unified, contextual messaging utilities intended for the 'tidyomics' ecosystem.

License GPL-3

License_use_gpl3 file LICENSE

Encoding UTF-8

Imports magrittr, tidyr, tibble, dplyr, purrr, pillar, S4Vectors, stringr, cli, vctrs, SummarizedExperiment, methods, rlang, fansi, pkgconfig

Roxygen list(markdown = TRUE)

RoxygenNote 7.3.3

Suggests knitr, rmarkdown, BiocStyle, airway, pasilla, testthat (>= 3.0.0)

Config/testthat/edition 3

VignetteBuilder knitr

URL <https://github.com/tidyomics/tidyprint>

BugReports <https://github.com/tidyomics/tidyprint/issues>

biocViews Software, Visualization, Infrastructure

Config/pak/sysreqs libicu-dev zlib1g-dev

Repository <https://bioc.r-universe.dev>

Date/Publication 2026-04-28 13:06:26 UTC

RemoteUrl <https://github.com/bioc/tidyprint>

RemoteRef HEAD

RemoteSha 9177293cacdfea376105d19dd0b4200e3c2dd542

Contents

tidyprint-package	2
as_tibble.SummarizedExperiment	3
demo_tidy_message	5
format_covariate_header	5
print.SummarizedExperiment	6
se_airway	7
tbl_format_header.SE_print_abstraction	8
tidy_message	8
tidy_print_toggle	9

Index **11**

tidyprint-package *tidyprint: Custom Print Methods for SummarizedExperiment*

Description

Provides customized print methods for [SummarizedExperiment](#) objects to enhance readability and usability within a tidy workflow. It offers consistent, tidyverse-aligned console displays, including alternative tibble abstractions for large genomic data to improve discoverability and interpretation. The package also includes unified, contextual messaging utilities intended for the 'tidyomics' ecosystem.

Value

This package provides S3 methods and functions that enhance the display and manipulation of [SummarizedExperiment](#) objects. The main functions return:

- `print.SummarizedExperiment`: Returns the input object invisibly
- `as_tibble.SummarizedExperiment`: Returns a tibble with combined assay and metadata
- `tidy_message`: Returns NULL invisibly (prints styled messages)

Key Features

- **Custom Print Methods**: Alternative display formats for [SummarizedExperiment](#) objects
- **Tidyverse Integration**: Seamless integration with tidyverse workflows
- **Scalable Visualization**: Adaptive displays for large genomic datasets
- **Standardized Messaging**: Consistent messaging utilities for the tidyomics ecosystem

Print Designs

The package provides multiple print designs for SummarizedExperiment objects:

- "SummarizedExperiment": Standard Bioconductor summary display
- "tidyprint_1": Compact tibble abstraction with adaptive separator rows
- "tidySummarizedExperiment": Legacy tidySummarizedExperiment style
- "plyxp": Compact top/bottom row display

Messaging System

The `tidy_message` function provides standardized, visually appealing messages for packages within the tidyomics ecosystem. It automatically detects the calling package to provide contextualized messaging.

Data Conversion

The `as_tibble` method converts SummarizedExperiment objects to tibble format, combining assay data with sample and feature metadata in a long format suitable for tidyverse workflows.

Author(s)

Chen Zhan <chen.zhan@adelaide.edu.au>

References

tidyprint: Custom Print Methods for SummarizedExperiment <https://github.com/tidyomics/tidyprint>

See Also

- [SummarizedExperiment](#) - The core Bioconductor class
- [as_tibble](#) - Base tibble conversion methods
- [print.SummarizedExperiment](#) - Custom print methods
- [tidy_message](#) - Messaging utilities

as_tibble.SummarizedExperiment

Convert SummarizedExperiment to tibble

Description

Converts a [SummarizedExperiment](#) object to a tibble format, combining assay data with sample and feature metadata in a long format suitable for tidyverse workflows.

Usage

```
## S3 method for class 'SummarizedExperiment'
as_tibble(
  x,
  ...,
  .name_repair = c("check_unique", "unique", "universal", "minimal"),
  rownames = pkgconfig::get_config("tibble::rownames", NULL)
)
```

Arguments

x	A SummarizedExperiment object to convert.
...	Additional arguments passed to internal conversion functions.
.name_repair	Treatment of problematic column names. See as_tibble for details.
rownames	Either NULL or a string giving the name of a column to use as rownames. See as_tibble for details.

Details

This method provides a bridge between Bioconductor's SummarizedExperiment objects and tidyverse data manipulation workflows. The conversion creates a long-format tibble where each row represents a feature-sample combination, making it suitable for filtering, grouping, and other tidyverse operations.

Value

A tibble containing the assay data combined with sample and feature metadata. The structure includes:

- Feature identifiers (from rownames or rowData)
- Sample identifiers (from colnames or colData)
- Assay values (one column per assay)
- Sample metadata (from colData)
- Feature metadata (from rowData)

See Also

[SummarizedExperiment](#), [as_tibble](#)

Examples

```
## Not run:
library(tidyprint)
data(se_airway)
as_tibble(se_airway)

## End(Not run)
```

demo_tidy_message	<i>Demonstrate tidy_message usage</i>
-------------------	---------------------------------------

Description

Runs through examples of the `tidy_message()` function for all supported message types ("info", "success", "warning", "danger"). Intended for demonstration and testing only.

Usage

```
demo_tidy_message()
```

Value

Invisibly returns NULL. Called for side effects (messages printed to the console).

Examples

```
demo_tidy_message()
```

format_covariate_header	<i>Format covariate header by distributing label across covariate columns</i>
-------------------------	---

Description

Format covariate header by distributing label across covariate columns

Usage

```
format_covariate_header(  
  separator_row,  
  printed_colnames,  
  covariate_names,  
  number_of_total_rows,  
  label = "COVARIATES"  
)
```

Arguments

<code>separator_row</code>	The separator row with column widths
<code>printed_colnames</code>	The printed column names
<code>covariate_names</code>	The names of covariate columns
<code>number_of_total_rows</code>	The total number of rows for spacing
<code>label</code>	The label to distribute (default: "COVARIATES")

Value

Formatted header string

```
print.SummarizedExperiment
```

Print method for SummarizedExperiment with tidyprint styles

Description

Provides a compact tibble-like display for a [SummarizedExperiment](#) object that preserves assay values and key covariates with a separator band when the table is truncated.

Usage

```
## S3 method for class 'SummarizedExperiment'
print(x, n = 10, ...)
```

Arguments

<code>x</code>	A <code>SummarizedExperiment</code> object to print.
<code>n</code>	Integer (default 10). Approximate number of rows to show in the display. When the total cells shown are fewer than <code>n</code> , the full table is printed and the separator row is suppressed.
<code>...</code>	Additional arguments passed to internal printers (currently unused).

Details

This method constructs a tibble abstraction for `SummarizedExperiment` data with columns: `.feature`, `.sample`, assay columns, a vertical separator `"|"`, followed by selected `colData` and `rowData` fields. When the output is truncated, an adaptive dash-only separator row is inserted after the first half block of rows. Additional indication of `colData` is provided as well.

By default, this method uses the standard `SummarizedExperiment` print format. To enable the tidy print format, use `tidy_print_on`. To check the current status, use `tidy_print_enabled`.

Value

x is returned *invisibly* after printing.

See Also

[SummarizedExperiment](#), [as_tibble](#)

Examples

```
## Not run:
library(tidyprint)
print(se_airway)                # compact tibble display
print(se_airway, n = 20)        # with custom row count

## End(Not run)
```

se_airway

Example SummarizedExperiment dataset

Description

Scaled RNA-seq counts for airway samples packaged as a [SummarizedExperiment](#).

Usage

```
data(se_airway)
```

Format

A SummarizedExperiment object with:

assays A matrix of scaled counts

colData Sample metadata (cell type, treatment, etc.)

Source

airway package

```
tbl_format_header.SE_print_abstraction
  Custom header for SE_print_abstraction
```

Description

Formats the header using pillar's formatting context so it aligns with the tibble body and adds a centered COVARIATES banner.

Usage

```
## S3 method for class 'SE_print_abstraction'
tbl_format_header(x, setup, ...)
```

Arguments

x	An object of class SE_print_abstraction.
setup	A pillar formatting setup (list) passed by pillar; contains things like the rendering width and the pre-rendered first body line (setup\$body).
...	Unused; present for S3 method consistency.

Value

A character vector of header lines to be displayed above the table. (pillar prints the returned character; this function has no side effects beyond returning those lines.)

```
tidy_message      Print a styled tidyprint message
```

Description

Prints a message to the console with a consistent tidyomics style. A package-specific prefix is automatically added based on the calling package name, followed by the message text. Different message types are displayed using appropriate styles from the **cli** package.

Usage

```
tidy_message(message, type = c("info", "success", "warning", "danger"))
```

Arguments

message	A character string containing the message to display.
type	The type of message to display. One of "info", "success", "warning", or "danger". Defaults to "info".

Value

Invisibly returns NULL. Called for side effects (printing a styled message to the console).

Examples

```
tidy_message("Loading data...", type = "info")
tidy_message("Data loaded successfully!", type = "success")
```

tidy_print_toggle	<i>Enable or disable tidy print for SummarizedExperiment objects</i>
-------------------	--

Description

These functions control whether SummarizedExperiment objects use the custom tidy print format or the standard Bioconductor print format. By default, standard print is used.

Usage

```
tidy_print_on(remember = FALSE)

tidy_print_off(remember = FALSE)

tidy_print_enabled()
```

Arguments

remember	Logical (default FALSE). If TRUE, saves the setting to a local cache file so it persists across R sessions. If FALSE, the setting only affects the current R session.
----------	---

Details

The tidy print format provides a compact tibble-like display that combines assay data with sample and feature metadata. The standard print format is the default Bioconductor display.

When tidy print is enabled, SummarizedExperiment objects will display using the custom tibble abstraction format. When disabled (default), they will use the standard SummarizedExperiment print method from the Bioconductor package.

The setting is stored in the global R options as `tidyprint.use_tidy_print`. When `remember = TRUE`, the setting is also saved to a cache file in the user's R configuration directory, which takes precedence over the option. The cache file location is determined by `tools::R_user_dir("tidyprint", "config")`.

Value

- `tidy_print_on()`: Returns TRUE invisibly after enabling tidy print
- `tidy_print_off()`: Returns FALSE invisibly after disabling tidy print
- `tidy_print_enabled()`: Returns a logical indicating whether tidy print is currently enabled

See Also

[print.SummarizedExperiment](#)

Examples

```
# Check current status
tidy_print_enabled()

# Enable tidy print (session only)
tidy_print_on()

# Enable tidy print and remember the setting (saved to cache)
tidy_print_on(remember = TRUE)

# Disable tidy print (use standard print)
tidy_print_off()

# Disable tidy print and remember the setting
tidy_print_off(remember = TRUE)
```

Index

* datasets

- se_airway, [7](#)

- as_tibble, [3](#), [4](#), [7](#)
- as_tibble.SummarizedExperiment, [3](#)

- demo_tidy_message, [5](#)

- format_covariate_header, [5](#)

- print.SummarizedExperiment, [3](#), [6](#), [10](#)

- se_airway, [7](#)
- SummarizedExperiment, [2-4](#), [6](#), [7](#)

- tbl_format_header.SE_print_abstraction,
[8](#)
- tidy_message, [3](#), [8](#)
- tidy_message(), [5](#)
- tidy_print_enabled, [6](#)
- tidy_print_enabled(tidy_print_toggle),
[9](#)
- tidy_print_off(tidy_print_toggle), [9](#)
- tidy_print_on, [6](#)
- tidy_print_on(tidy_print_toggle), [9](#)
- tidy_print_toggle, [9](#)
- tidyprint(tidyprint-package), [2](#)
- tidyprint-package, [2](#)