Package: qcmetrics (via r-universe)

	June 30, 2024
Type Pac	kage
Title AF	ramework for Quality Control
Version 1	1.43.0
con indi	on The package provides a framework for generic quality trol of data. It permits to create, manage and visualise vidual or sets of quality control metrics and generate lity control reports in various formats.
Depends	R (>= 3.3)
Imports	Biobase, methods, knitr, tools, xtable, pander, S4Vectors
	affy, MSnbase, ggplot2, lattice, mzR, BiocStyle, rmarkdown, kdown
URL htt	p://lgatto.github.io/qcmetrics/articles/qcmetrics.html
BugRepo	<pre>rts https://github.com/lgatto/qcmetrics/issues</pre>
	s ImmunoOncology, Software, QualityControl, Proteomics, croarray, MassSpectrometry, Visualization, ReportWriting
VignetteI	Builder knitr
Roxygenl	Note 6.1.0
Repositor	ry https://bioc.r-universe.dev
RemoteU	rl https://github.com/bioc/qcmetrics
RemoteR	ef HEAD
RemoteS	ha efceb1344d3015ade5435cbaa9ba00c64161a272
Conte	nts
1	n15qc psm Qc2Tex QcMetadata-class QcMetric-class QcMetrics-class QcMetrics-class

n15qc

Index 11

n15qc N15 labelling QC report	
25 10001	

Description

A simple wrapper for the QC of 15N labelling. The respective QC items are the distribution of PSM incorporation rates, distribution of log2 fold-changes and number of identified features. See the vignette for details.

Usage

```
n15qc(object,
fcol = c("Protein_Accession", "Peptide_Sequence", "Number_Of_Unique_Peptides", "Variable_Modificat
inctr = 97.5, lfctr = c(-0.5, 0.5), type, reportname)
```

Arguments

object	An MSnSet to be quality controlled.
fcol	The name of the feature variables for the protein identifiers (accession numbers for example), the peptide sequences, the number of unique peptides for each identified protein, the variable modifications identified on the peptides and the N15 incorporation rate. These must be provided in that order. Defaults are Protein_Accession, Peptide_Sequence, Number_Of_Unique_Peptides, Variable_Modifications, and inc.
inctr	The minimum level of median incorporation rate to set the QC item status to TRUE. Default is 97.5.
lfctr	The range of accepted median PSM $\log 2$ fold-change for the QC item status to be set to TRUE. Default is $c(-0.5, 0.5)$.
type	The type of report to be saved. If missing (default), no report is generated. See qcReport for details.
reportname	The name of the report, in case a type is defined. If missing (default), the report will be names n15qcreport followed by the generation data and time.

Value

Invisibly returns the resulting QcMetrics instance.

Author(s)

Laurent Gatto

psm 3

psm

15N example data

Description

An example data for 15N metabolic labelling, distributed as an MSnSet to illustrate 15N QC.

Usage

```
data(n15psm)
```

References

See MSnSet and the MSnbase-demo vignette in the MSnbase package.

Examples

```
library("MSnbase")
data(n15psm)
psm
```

Qc2Tex

'QcMetric' sectioning functions

Description

These functions convert the ith QcMetric instance of the QcMetrics object into a section of the adequate format, i.e. TeX of R markdown.

Usage

```
Qc2Tex(object, i)
Qc2Tex2(object, i)
Qc2Tex3(object, i)
Qc2Rmd(object, i)
```

Arguments

object An instance of class QcMetrics with at least one QC item.

i A numeric of length 1 indicating the index of the item to be conve

A numeric of length 1 indicating the index of the item to be converted into text section.

Value

A character representing the QC item section.

4 QcMetadata-class

Author(s)

Laurent Gatto

See Also

qcReport and the vignette.

QcMetadata-class

The "QcMetadata" class

Description

The QcMetadata class is a simple interface to metadata. The objects can be displayed with show for a summary and print for the content.

Objects from the Class

Objects can be created by calls of the form QcMetadata(...).

Slots

metadata: Object of class "list" that stores the metadata variables. The list must be named. NA and empty characters are not allowed.

Methods

```
[ signature(x = "QcMetadata"): subsets x as a new QcMetadata instance.
```

[[signature(x = "QcMetadata"): extracts a single element of x.

metadata signature(x = "QcMetadata"): return the object's metadata list. Also available as mdata.

metadata<- signature(x = "QcMetadata", value = "list"): sets the objects metadata. Also
 available as mdata.</pre>

length signature(x = "QcMetadata"): returns the number of metadata variables.

names signature(x = "QcMetadata"): returns the names of the metadata variables.

Author(s)

Laurent Gatto

QcMetric-class 5

Examples

QcMetric-class

The "QcMetric" class for QC items

Description

Data structure for individual QC metrics items.

Objects from the Class

Objects can be created using the constructor QcMetric(...), where slots are assigned individually. See example below.

Slots

name: Object of class "character" of length 1 naming the object.

description: Object of class "character" of arbitrary length describing the qc metric in more details.

qcdata: Object of class "environment" that stores the actual data.

plot: Object of class "function" to graphically represent the data and infer quality status.

show: Object of class "function" to print a short textual representation of the object. A reasonable default value is provided.

status: Object of class "logical" that indicates weather the data passes (TRUE) or fails (FALSE) the metric or has not yet been evaluated.

Methods

```
name signature(object = "QcMetric"): retrieves the name of the object.
name<- signature(object = "QcMetric", value = "character"): set the name of the object.
description signature(object = "QcMetric"): retrieves the description of the object.
description<- signature(object = "QcMetric", value = "character"): set the description of the object.
status signature(object = "QcMetric"): retrieves the status of the object.</pre>
```

status<- signature(object = "QcMetric", value = "logical"): sets the status of the objects.

6 QcMetric-class

```
qcdata signature(object = "QcMetric", x = "missing"): lists all the data objects that are
associated with the objects.
```

qcdata signature(object = "QcMetric", x = "character"): retrieves the variable x for the object.

qcdata<- signature(object = "QcMetric", var): creates or overwrites (with a message) the data variable var by assigning the RHS value. If var is missing and the RHS expression is an environment, then qcdata is reset with all the variables in value.

qcenv signature(object = "QcMetric"): return the environment that stores the QC data.

qcenv<- signature(object = "QcMetric"): Set all variable in the RHS environment as qcdata variables. Equivalent to qcdata(object) <- x where x is an environment.

show signature(object = "QcMetric"): shows a textual summary of object. The default show implementation is available as the qcshow{object, qcdata} function. The second argument is a logical (default is TRUE) that specifies whether qcdata(object) should be displayed.

show<- signature(object = "QcMetric", value = "function"): sets a custom show method
for object.</pre>

plot signature(x = "QcMetric", y = "missing"): plots the object using the provide show method.

plot<- signature(object = "QcMetric", value = "function"): sets a custom plot method for object.

qcReport signature(x = "QcMetric", ...): to generate quality reports. See qcReport for details.

Author(s)

Laurent Gatto

See Also

The QcMetrics class to bundle a set of QcMetric instances.

Examples

```
(qc <- QcMetric())</pre>
qcdata(qc)
try(qcdata(qc, "x"))
x \leftarrow rnorm(10)
qcdata(qc, "qc1") <- x
qcdata(qc, "qc2") <- 1:10</pre>
qcdata(qc)
all.equal(qcdata(qc, "qc1"), x)
all.equal(qcdata(qc, "qc2"), 1:10)
name(qc) <- "My test QcMetric"</pre>
description(qc) <- "This qc metric describes bla bla bla, indicating possible issues in the third step of protocol b
status(qc) <- FALSE
qc
## or
e <- new.env()
e$qc1 <- rnorm(100)
```

QcMetrics-class 7

```
e$qc2 <- 1:100
qcenv(qc) <- e
length(qcdata(qc, "qc1"))
head(qcdata(qc, "qc2"))
show(qc)
show(qc) <- function(object) cat("Updated show method\n")</pre>
show(qc)
show(qc) <- qcshow
qc
plot(qc)
plot(qc) <-
    function(object, ...)
        plot(qcdata(object, "qc2"),
             qcdata(object, "qc1"),
             xlab = "qc1",
             ylab = "qc2",
              ...)
plot(qc)
plot(qc, col = "red", pch = 19)
```

QcMetrics-class

The "QcMetrics" class for collections of QC items

Description

Data structure for storing lists of QcMetric items.

Objects from the Class

Objects can be created using the constructor QcMetrics(...), where slots are assigned individually. See example below.

In a standardised quality control pipeline, the QcMetrics and QcMetric object are not generated manually. Their creation is delegated to a wrapper function that reads a specific type of files, parses the data, produces the individual QcMetric instances and, eventually, the QcMetric object. See the package vignette for details and examples.

Slots

metadata: Object of class QcMetadata storing the metadata of the object. This list would typically contain the input file the data was read from, the date the object was generated, ... or fully fledged *minimum information* descriptions (see MIAXE), when available.

qcdata: Object of class "list" storing all the individual QcMetric instances.

8 QcMetrics-class

Methods

```
[ signature(x = "QcMetrics"): subsets x as a new QcMetrics instance.
[[ signature(x = "QcMetrics"): extracts a single QcMetric instance.
length signature(x = "QcMetrics"): returns the number of QcMetric instances populate x.
metadata signature(x = "QcMetrics"): return the object's metadata list. Also available as
metadata<- signature(x = "QcMetrics", value = "list"): sets the objects metadata. Also</pre>
     available as mdata.
metadata<- signature(x = "QcMetric", value = "QcMetadata"): sets the objects metadata. Also
     available as mdata.
name signature(object = "QcMetrics"): returns a character vector of length length(object)
     with the names of the QcMetric instances.
qcdata signature(object = "QcMetrics", x = "missing"): returns a list of all QcMetric in-
     stances.
qcdata<- signature(object = "QcMetrics", value = "list"): sets the qcdata of object.</pre>
show signature(object = "QcMetrics"): prints a short textual description of object.
status signature(object = "QcMetrics"): returns a vector of quality statuses (logicals).
status<- signature(object = "QcMetrics", value = "logical"): sets the quality statuses. Length
     of value and object must be identical.
as signature(object = "QcMetrics", "data.frame"): coerces object as a length(object)
     by 2 data frame with the respective QcMetric instances names and statuses.
qcReport signature(object = "QcMetrics"): ...
```

Author(s)

Laurent Gatto

See Also

The OcMetric class for individual QC items.

Examples

qcReport-methods 9

qcReport-methods

Automatic report generation

Description

The qcReport method generates report in various formats taking a QcMetrics instance as input. Each individual quality control item produces a section with description of the item and a assessment figure.

Details

The reporting functions take a QcMetrics instance as input, generate the source of the report and compile it into the final format that are currently available are reporting_pdf, reporting_tex, reporting_rmd and reporting_html. See qcto for details about the sectioning functions, that convert individual QcMetric objects into adequate report sections.

The package vignette documents the report generation in more details and describes possibilities for customisation.

Methods

signature(object = "QcMetrics", reportname = "character", type = "character", author = "character", title
generates a report for the QcMetrics object. The report will be named according the reportname
(default is qcreprt)and type, the latter defining the output format and the extension. Possible types are pdf (default), "tex", "Rmd", "html" (all generated using the package knitr). A
custom title can be provided; default is "Quality control report generated with qcmetrics".
If no author is provided, the default value (Sys.getenv("USER")) is used. The addition of
a table of contents (default is FALSE), a metadata section, a summary section and the session
information can be controlled with the toc, metadata, summary and sessioninformation
arguments. The metadata section is added to the report when present and the other have TRUE
as default.

10 qcReport-methods

Intermediate files are deleted, unless clean is set to FALSE and verbose output can be turned on by setting quiet to FALSE.

The reporter and qcto arguments are used to convert QcMetric and QcMetrics objects into report source. See Details and the package vignette for details.

Addition parameters can be passed to inner functions. For the pdf report, passed to texi2pdf; for html, passed to markdown::markdownToHTML.

The method invisibly returns the name of the report that was generated.

Examples

```
example(QcMetrics)
show(qcm)

destdir <- tempdir()
(report <- file.path(destdir, "testQCReport"))

## Not run:
    ## pdf report
    qcReport(qcm, reportname = report)
    ## use pdflatex to generate the pdf file
    qcReport(qcm, reportname = report, texi2dvi = "pdflatex")

## End(Not run)

## default html report
html <- qcReport(qcm, reportname = report, type = "html")
html
if (interactive())
    browseURL(html)</pre>
```

Index

* classes	mdata<- (QcMetrics-class), 7	
QcMetadata-class, 4	<pre>mdata<-,QcMetadata,list-method</pre>	
QcMetric-class, 5	(QcMetadata-class), 4	
QcMetrics-class, 7	<pre>mdata<-,QcMetrics,list-method</pre>	
* datasets	(QcMetrics-class), 7	
psm, 3	mdata<-,QcMetrics,QcMetadata-method	
* methods	(QcMetrics-class), 7	
qcReport-methods, 9	metadata(QcMetrics-class),7	
* report	metadata,QcMetadata-method	
qcReport-methods, 9	(QcMetadata-class), 4	
[,QcMetadata-method(QcMetadata-class),	metadata,QcMetrics-method	
4	(QcMetrics-class), 7	
[,QcMetrics-method(QcMetrics-class),7	<pre>metadata<- (QcMetrics-class), 7</pre>	
[[,QcMetadata-method	<pre>metadata<-,QcMetadata,list-method</pre>	
(QcMetadata-class), 4	(QcMetadata-class), 4	
[[,QcMetrics-method(QcMetrics-class),7	<pre>metadata<-,QcMetrics,list-method</pre>	
	(QcMetrics-class), 7	
as.data.frame.QcMetrics	metadata<-,QcMetrics,QcMetadata-method	
(QcMetrics-class), 7	(QcMetrics-class), 7	
	MIAxE, 7	
class:QcMetadata (QcMetadata-class), 4		
class:QcMetric (QcMetric-class), 5	n15psm (psm), 3	
class:QcMetrics(QcMetrics-class),7	n15qc, 2	
description,QcMetric-method	name (QcMetrics-class), 7	
(QcMetric-class), 5	name, QcMetric-method (QcMetric-class), 5	
description<-,QcMetric,character-method	name,QcMetrics-method	
(QcMetric-class), 5	(QcMetrics-class), 7	
(QCHeth IC Class), 3	<pre>name<- (QcMetrics-class), 7</pre>	
example_reports (qcReport-methods), 9	name<-,QcMetric,character-method	
2 1 P 2 2 3 P 2 1 2 (4 2 3 P 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2	(QcMetric-class), 5	
length,QcMetadata-method	name<-,QcMetrics,character-method	
(QcMetadata-class), 4	(QcMetrics-class), 7	
length,QcMetrics-method	names,QcMetadata-method	
(QcMetrics-class),7	(QcMetadata-class), 4	
	names<-,QcMetadata,character-method	
mdata(QcMetrics-class),7	(QcMetadata-class), 4	
mdata,QcMetadata-method		
(QcMetadata-class), 4	plot,QcMetric,missing-method	
mdata,QcMetrics-method	(QcMetric-class), 5	
(QcMetrics-class), 7	<pre>plot<- (QcMetric-class), 5</pre>	

12 INDEX

plot<-,QcMetric,function-method	<pre>reporting_rmd (qcReport-methods), 9</pre>
(QcMetric-class), 5	reporting_tex (qcReport-methods), 9
print,QcMetadata-method	
(QcMetadata-class), 4	show,QcMetadata-method
psm, 3	(QcMetadata-class), 4
,	show, QcMetric-method (QcMetric-class), 5
Qc2Rmd (Qc2Tex), 3	show,QcMetrics-method
Qc2Tex, 3	(QcMetrics-class), 7
Qc2Tex2 (Qc2Tex), 3	<pre>show<- (QcMetric-class), 5</pre>
Qc2Tex3 (Qc2Tex), 3	<pre>show<-,QcMetric,function-method</pre>
qcdata(QcMetric-class), 5	(QcMetric-class), 5
qcdata, QcMetric, character-method	status (QcMetric-class), 5
(QcMetric-class), 5	status,QcMetric-method
qcdata,QcMetric,missing-method	(QcMetric-class), 5
(QcMetric-class), 5	status,QcMetrics-method
qcdata,QcMetrics,ANY-method	(QcMetrics-class), 7
(QcMetrics-class), 7	status<- (QcMetric-class), 5
qcdata<- (QcMetric-class), 5	status<-,QcMetric,logical-method
qcdata<-,QcMetric,ANY-method	(QcMetric-class), 5
(QcMetric-class), 5	status<-,QcMetrics,logical-method
qcdata<-,QcMetrics,list-method	(QcMetrics-class), 7
(QcMetrics-class), 7	(2
qcenv (QcMetric-class), 5	
qcenv,QcMetric-method (QcMetric-class),	
5	
qcenv<- (QcMetric-class), 5	
qcenv<-,QcMetric,environment-method	
(QcMetric-class), 5	
OcMetadata, 7	
QcMetadata, / QcMetadata (QcMetadata-class), 4	
QcMetadata-class, 4	
QcMetric, 7, 8	
QcMetric (QcMetric-class), 5	
QcMetric-class, 5	
QcMetrics, 6, 9	
QcMetrics (QcMetrics-class), 7	
QcMetrics-class, 7	
qcReport, 2, 4, 6	
qcReport (qcReport-methods), 9	
qcReport,QcMetrics-method	
(qcReport-methods), 9	
qcReport-methods, 9	
qcshow (QcMetric-class), 5	
qcto, 9	
qcto (Qc2Tex), 3	
reporting (qcReport-methods), 9	
<pre>reporting_html (qcReport-methods), 9</pre>	
reporting pdf (acReport-methods), 9	